

# Andrasko Staining Grid

Percentage of Average Corneal Staining Area at 2 Hours

		Branded Solutions										Private Label Solutions			
		Unisol <sup>1</sup> 4 Saline	Clear Care <sup>4</sup>	OPTI-FREE <sup>®</sup> EXPRESS <sup>1</sup>	OPTI-FREE <sup>®</sup> RepleniSH <sup>1</sup>	OPTI-FREE <sup>®</sup> PureMoist <sup>1</sup>	Biotrue <sup>3</sup>	Renu Fresh <sup>3</sup>	Renu Sensitive <sup>3</sup>	Complete MPS <sup>2</sup>	Aquify <sup>4</sup>	Walmart MPS (Renu M+)	Target MPS (Renu M+)	CVS MPS (Renu M+)	Walgreen MPS (Renu M+)
Hydrogel	Acuvue <sup>5</sup> 2	1%	1%	2%	5%	1%	1%	1%	1%	1%	1%	1%	1%	1%	1%
	Proclear <sup>6</sup>	1%	1%	1%	2%	1%	28%	57%	23%	6%	12%	61%	54%	53%	42%
	Soflens <sup>3</sup> 66	1%	1%	1%	1%	1%	52%	73%	32%	17%	8%	66%	62%	63%	56%
Silicone-Hydrogels	Acuvue Advance <sup>5</sup>	1%	1%	1%	1%	1%	9%	13%	4%	12%	2%	16%	13%	12%	12%
	Acuvue Oasys <sup>5</sup>	2%	1%	3%	5%	2%	1%	9%	5%	4%	3%	12%	8%	13%	10%
	Biofinity <sup>6</sup>	2%	2%	3%	2%	1%	17%	4%	2%	2%	2%	4%	3%	3%	2%
	Purevision <sup>3</sup>	2%	1%	4%	7%	3%	46%	73%	43%	15%	21%	71%	76%	No Testing Planned	No Testing Planned
	O2 Optix <sup>4</sup>	2%	1%	2%	5%	1%	21%	24%	7%	3%	3%	41%	28%	28%	24%
	Night & Day <sup>4</sup>	2%	1%	2%	3%	1%	17%	24%	11%	1%	3%	36%	24%	26%	22%
Updated: August 19, 2011		H <sub>2</sub> O <sub>2</sub>	POLYQUAD/ALDOX				PHMB/ Polyquaternium		BIGUANIDES (PHMB)						

Staining Zone Color Codes

■ under 10%    ■ 10% to 20%    ■ over 20%

Trademarks: 1 Alcon; 2 AMO; 3 Bausch & Lomb; 4 Novartis; 5 Johnson & Johnson; 6 CooperVision

Data Courtesy of: [www.StainingGrid.com](http://www.StainingGrid.com)